

# Edwards General Surgery Catheters

A Spectrum of Solutions for Identification and Removal of Gallstones

Edwards Lifesciences offers a range of products for the identification and removal of gallstones, including cholangiography and laparoscopic cholangiography catheters, cholangiography kits, Fogarty biliary balloon probe, and Fogarty irrigation catheters.

## Cholangiography Catheter

- Allows introduction of contrast media into cystic and common bile ducts
- Clear vinyl body makes it easy to detect air bubbles, which could be mistaken for stones
- Three-way luer lock hub enables the introduction of irrigation fluid and contrast media
- Black marking at the distal tip extends 1.5 cm proximally and facilitates placement
- Beveled tip facilitates introduction into common duct and reduces the risk of accidental withdrawal

## Laparoscopic Cholangiography Catheter

- Cholangiography catheter also available in 80 cm length to meet the needs for laparoscopic procedures



# Edwards General Surgery Catheters

## A Spectrum of Solutions for Identification and Removal of Gallstones

### Cholangiography Kit

- Specialized instruments needed for operative cholangiography and “mini” cholecystectomy procedures provided in a single-pack sterile, disposable kit:
  - 6 mm spring clip and spring clip applicator
  - 4F cholangiography catheter
  - two 20 cc syringes
- Unique spring clip port allows for maximum injection of contrast media without accidental withdrawal
- Spring clip is maneuverable in 30° increments up to 180° and allows the necessary angle for placement

### Fogarty Biliary Balloon Probe

- Balloon catheter designed to atraumatically remove hepatic calculi
- Optimally thick balloon resists puncture by the calculi
- Radiopaque body is stiff enough to manipulate, yet flexible enough to minimize damage to the duct wall

### Fogarty Irrigation Catheter

- 10F catheter designed to facilitate flushing contrast media from the ductal system

SPECIFICATIONS			
Model Number	Catheter French Size	Length	Maximum Liquid Capacity
<b>Cholangiography Catheter</b>			
420404F	4F	40 cm	
420406F	6F	40 cm	
<b>Laparoscopic Cholangiography Catheter</b>			
420804F	4F	80 cm	
<b>Cholangiography Kit</b>			
CV8000 (includes catheter, clip, applicator, 2 syringes)	4F	40 cm	
<b>Fogarty Biliary Balloon Probe</b>			
410235F	5F	23 cm	0.4 ml
410405F	5F	40 cm	0.4 ml
410236F	6F	23 cm	1.25 ml
<b>Fogarty Irrigation Catheter</b>			
220804F	4F	80 cm	
220806F	6F	80 cm	
2202310F	10F	23 cm	

*When you are seeking time-tested solutions for identification and removal of gallstones, reach for Edwards' general surgery catheters.*

For professional use. CAUTION: Federal (United States) law restricts this device to sale by or on the order of a physician. See instructions for use for full prescribing information, including indications, contraindications, warnings, precautions and adverse events.

Edwards is a trademark of Edwards Lifesciences Corporation. Edwards Lifesciences, the stylized E logo and Fogarty are trademarks of Edwards Lifesciences Corporation and are registered in the United States Patent and Trademark Office. © 2010 Edwards Lifesciences Corporation. All rights reserved. AR05677

Edwards Lifesciences

Irvine, USA | Nyon, Switzerland | Tokyo, Japan | Singapore, Singapore | Sao Paulo, Brazil

[www.edwards.com](http://www.edwards.com)

